



# RNAseq applications




# Our TWIST portfolio covers RNAseq solutions

- Targeted sequencing (RNA exome or small panel)
- Whole Transcriptome Sequencing


## TWIST TARGETED SEQUENCING WORKFLOW




## TWIST WHOLE TRANSCRIPTOME SEQUENCING WORKFLOW



# Twist RNAseq portfolio resolution



# We will find the right RNAseq approach for you



- **Twist Whole Transcriptome**
  - Deep sequence total mRNA to discover new transcripts and analyze the transcriptome
- **Twist RNA Exome**
  - Capture the protein coding sequencing for differential expression and novel fusion detection
- **Twist Alliance CeGaT Fusion Panel or Custom Panels**
  - Build custom capture panels for RNA libraries for targeted RNA analysis

